

**Product Name :**  
Direct Shear (Motorized 1 Speed)

**Product Code :**  
NLAB-TESTINGLABM109047



**Description :**

Direct Shear (Motorized 1 Speed)

**Technical Specification :**

Same as above but electrically operated to give a single constant rate of stain of 1.25mm/min. Suitable for 230V A.C. Weight 175kg. approx.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003