

**Product Name :**  
Disc Type Wire Gauge

**Product Code :**  
NLAB-LBEMALAB230030



**Description :**

Disc Type Wire Gauge

**Technical Specification :**

The Wire Gauge can measure from 1 to 36 and hence is most versatile in application

Made of Hardened and Tempered Steel adopted from Japan, which ensures no wear and tear of the gauge during measurements

This gauge is a boon for Wire manufactures, Plate manufacturers, Fabricators, Traders and Jewellers

The preciseness of the dimensions has made wire gauge a very popular Measuring Tool among Jewellers

Wire Gauge is a very handy and day to day usage precise tool to measure thickness of Metal Plate or Dia of a Wire

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

