

Product Name :
Displacement Measurement (LVDT)

Product Code :
NLAB-ENGINEERINGLB52001



Description :

Displacement Measurement (LVDT)

Technical Specification :

DESCRIPTION

Here simple demonstration setup is provided to experiment with LVDT for measurement of linear displacement. To vary the actual displacement an annual micrometer driven mechanism will be provided. Whereas for reference comparison same micrometer will be used.

SPECIFICATIONS

1. LVDT sensor with necessary signal conditioning circuit.
2. Micrometer comparison scale.
3. Manual micrometer driven mechanism for displacement variation.
4. Digital Indicator: -

Least Count: - 0.1 mm

Digit: - 3 1/2

5. Range – 0 to 20 mm.

6. Powder coated enclosure

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003