

Product Name :
Dissolution Rate Test Equipment

Product Code :
NLAB-LABORATORYEM65019



Description :

Naugra Lab Dissolution Rate Test Equipment

Technical Specification :

Dissolution Rate Test Equipment:-

This equipment is used to determine dissolution rate of tablets and complies with US-Pharmacopoeia-XVII standard with water constantly moving and holding temp.

Constant at $37^{\circ}\text{C} \pm 0.5^{\circ}\text{C}$ throughout the test.

A 1000 ml capacity cylindrical vessel with flanges near the top is provided to accept 4 parts cover.

Motor speed is controlled by a speed controller within a range of 25 to 150 RPM.

A stainless steel basket is provided to hold samples.

The machine works on 220V AC 50Hz single phase.

Single test

Combined unit for 3 tests.

Combined unit for 6 tests.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003