

Product Name :
Dissolved Oxygen Meters

Product Code :
NLAB-LABORATORYEM18002



Description :

Dissolved Oxygen Meters

Technical Specification :

Highly Stable and Accurate

DO measurements

3½ Digit LCD Display

ABS Plastic Body

Rugged Portability

Designed in ABS Plastic to provide the rugged portability matched in fields with the accuracy demanded in the laboratory.

A handy carrying case holds every thing needed in the field operation.

The instrument measures DO of any solution in the fields, laboratory or plant.

Range: 0 - 20 ppm

Resolution: 0.1 ppm

Accuracy: ± 0.1 ppm

Temp. Comp: 0 - 50 °C (Manual)

Calibration: Manual

Sensor: Gold/Silver Amperometric DO Probe

Display: 3½ Digit LCD

Power: 9 V DC Battery

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003