

**Product Name :**  
Dissolved Oxygen Meters

**Product Code :**  
NLAB-LABORATORYEM18002



**Description :**

Dissolved Oxygen Meters

**Technical Specification :**

Highly Stable and Accurate

DO measurements

3½ Digit LCD Display

ABS Plastic Body

Rugged Portability

Designed in ABS Plastic to provide the rugged portability matched in fields with the accuracy demanded in the laboratory.

A handy carrying case holds every thing needed in the field operation.

The instrument measures DO of any solution in the fields, laboratory or plant.

Range: 0 - 20 ppm

Resolution: 0.1 ppm

Accuracy:  $\pm 0.1$  ppm

Temp. Comp: 0 - 50 °C (Manual)

Calibration: Manual

Sensor: Gold/Silver Amperometric DO Probe

Display: 3½ Digit LCD

Power: 9 V DC Battery

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003