

Product Name :
Distance Sensor

Product Code :
NLAB-VOCATIONALTVET28002



Description :

Distance Sensor

Technical Specification :

Distance measuring range: 10 cm to 80 cm (4" to 32")

Output type: analog voltage

Output voltage differential over distance range: 1.9 V

Operating voltage: 4.5 V to 5.5 V

Average current consumption: 30 mA

Weight: 3.5 g (0.12 oz)

Response time: 38 ± 10 ms

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com