

**Product Name :**  
Distillation Column Simulator

**Product Code :**  
NLAB-ENGINEERINGLB36022



**Description :**

Distillation Column Simulator

**Technical Specification :**

Centrifugal pump

Feeding tank with 100 litres capacity.

Variable area flow meter with the glass tube

Centrifugal fan.

Electrical control panel in main cabinet.

Plexiglas distillation column diameter 600 mm with 2 sieve trays,

Emergency stop button.

A bubble cap tray and three pressure tappings.

2 differential manometers.

We are well-known manufacturers, OEM suppliers of Distillation Column Simulator for Chemical Engineering Lab Equipments. Contact us for high quality Distillation Column Simulator for Chemical Engineering Lab Equipments for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003