

Product Name :
Distillation Stills Barn stead Type

Product Code :
NLAB-EDUSCHOOLEQ35002



Description :

Distillation Stills Barn stead Type

Technical Specification :

- Barnstead Type Distillation Stills have table model, double walled construction with inner and outer chambers made of stainless steel sheet.
- Gap filled with glass wool for temp. Insulation.
- Supplied with plug-cord to work on **220V AC 50Hz** single phase.

Capacities available:

10 ltrs Per Hour
15 ltrs Per Hour
20 ltrs Per Hour
25 ltrs Per Hour

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003