

**Product Name :**  
Distillation Stills Manasty Type Wall Hanging Tools

**Product Code :**  
NLAB-EDUSCHOOLEQ35001



**Description :**

Distillation Stills Manasty Type Wall Hanging

**Technical Specification :**

- Unit made of all stainless steel sheet for making pyrogen free distilled water.
- Kettle element with cut-out system has been employed.
- Clamps have been provided for wall hanging.
- Supplied with plug-cord to work on **220V AC 50Hz** single phase.

**Capacities available:**

2 ltrs Per Hour  
4 ltrs Per Hour  
5 ltrs Per Hour  
6 ltrs Per Hour

8 ltrs Per Hour

16 ltrs Per Hour

20 ltrs Per Hour

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003