


Product Name : Donut Maker	Product Code : NLAB-VOCATIONALTVET27010
--------------------------------------	---



Description :

Donut Maker

Technical Specification :

Bake 6 light and fluffy doughnuts in about 4 minutes

Nonstick baking plates for ideal results and easy clean up

Easy to Use

Material: Stainless Steel

Voltage: 210 V

Frequency: 50 Hz

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003