

**Product Name :**  
Double Gel Electrophoresis Apparatus

**Product Code :**  
NLAB-EDUSCHOOLEQ21002



**Description :**

Double Gel Electrophoresis Apparatus

**Technical Specification :**

Run two gels simultaneously with the leak-free, the one-piece moulded chamber is constructed of durable, clear acrylic and has pure platinum electrodes for superior conductivity and longevity.

The system includes two casting gel trays and six 12-well forming combs and patch cords.

Specially designed casting trays hold melted agarose without the use of tape or snap-on ends.

Built-in safety features prevent the student from using the chamber unless the lid is locked in place.

Wells can be formed at the end or the centre of the gel.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003