

Product Name :
Double Gel Electrophoresis Apparatus

Product Code :
NLAB-EDUSCHOOLEQ21002



Description :

Double Gel Electrophoresis Apparatus

Technical Specification :

Run two gels simultaneously with the leak-free, the one-piece moulded chamber is constructed of durable, clear acrylic and has pure platinum electrodes for superior conductivity and longevity.

The system includes two casting gel trays and six 12-well forming combs and patch cords.

Specially designed casting trays hold melted agarose without the use of tape or snap-on ends.

Built-in safety features prevent the student from using the chamber unless the lid is locked in place.

Wells can be formed at the end or the centre of the gel.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003