

Product Name :
Double Juice Dispenser with Compressor

Product Code :
NLAB-VOCATIONALTVET27065

Description :

Double Juice Dispenser with Compressor

Technical Specification :

Machine Type: Automatic

Capacity: 18 Liters

Number of Flavours: 4

Body Material: Stainless Steel

Features:

Optimum performance

Sturdy body

Easy to clean

Quality tested

Highly efficient

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003