

Product Name :
Double Needle Lockstitch Machine

Product Code :
NLAB-VOCATIONALTVET22002



Description :

Double Needle Lockstitch Machine for Textile Vocational Technology Lab

Technical Specification :

Application: Light & Medium Material.

No of Needles: 2

Needle Gauge: 3/8"

Lubrication: Auto Lubrication.

Motor: 250W Clutch type with the reversible plug-in mechanism.

Maximum Stitch length: 0-4mm approx.

Presser Foot Lifting: 7-13mm approx.

Maximum Sewing Speed: 4000 SMP approx.

Stand & Table: Standard Table with Adjustable height Stand.

Standard Accessories & Video demo CD supplied with each machine

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003