

**Product Name :**  
Double Wall Autoclave With Vacuum

**Product Code :**  
NLAB-LABORATORYEM47012



**Description :**

Double Wall Autoclave With Vacuum

**Technical Specification :**

Double Wall Autoclave With Vacuum -

**Features:-**

Top end model for sterilization of all critical application with exact reproducibility

Unique double wall design

Fully automatic operation

µp based PID controller with dual display for Temp. and Time

Reduced loading height

Heavy duty jacket

Pre & Post Vacuum available

**Specifications:-**



500mm

Power  
Consumption

Operating /  
Pressure

Operating  
Temperature

230V

Digital

230kg

Weight

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003