

Product Name : Drager Evita Ventilators	Product Code : NLABMEDICAL-M3390011
Description : Drager Evita Ventilators India, Drager Evita Ventilators manufacturer India, suppliers, Medical Drager Evita Ventilators India, Hospital Drager Evita Ventilators India	
Technical Specification : Airway pressure: Peak pressure, plateau pressure, mean pressure, PEEP, min. pressure (0 to 99 mbar/cmH2O). Minute volume (MV), (BTPS): MV, MVspont, MVleak (0 to 99 L/min). Tidal volume (VT), (BTPS): Inspired VT, expired VT, VTASB (VTPS) (0 to 3999 mL). Airway pressure: High/Low. Expired minute volume: High/Low. Tidal volume: High. Apnea alarm time: 5 to 60 s. Spontaneous breath frequency: High . Max. flow for pressure support and spontaneous breathing: 180 L/min (adult), 60 L/min (pediatric). Mains power connection: 100 to 240 V, 50/60 Hz, 10 to 30 V DC (optional). Power consumption: Approx. 125 W. Gas supply operating pressure: O2, air: 2.7 to 6 bar/39 to 87 PSI. Dimensions ventilator (WxHxD): 530x290x450 cm (20.9x11.4x17.7 inches) (without trolley). Weight basic unit: Approx. 27 kg (60 lbs.). Diagonal screen size: 6.5" TFT LCD color display.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003