

Product Name :
Drager Oxylog 1000 Ventilator

Product Code :
NLABMEDICAL-M3390009

Description :

Drager Oxylog 1000 Ventilator

Technical Specification :

Oxylog 1000 – a time-cycled, volume controlled and pressure limited emergency ventilator for the controlled ventilation of patients who require a minute volume of at least 3 L/min.

Dimensions (W x H x D) 215 x 90 x 215 mm / 8.5 x 3.5 x 8.5 inches (excl. handle).

Weight 3.15 kg / 7.3 lbs.

Drive Gas: Medical grade O2 or in exceptional cases compressed air.

Supply pressure 2.7 to 6 bar / 40 to 88 psi at 60 L/min.

Performance Data:-Ventilation frequency, smoothly variable 4 to 54 1/min.

Minute volume, smoothly variable 3 to 20 L/min.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003