

<b>Product Name :</b> Drager Savina	<b>Product Code :</b> NLABMEDICAL-M3390010
--	---

--

**Description :**

Drager Savina

**Technical Specification :**

Ventilation modes - IPPV (CMV), IPPV Assist (CMV Assist).  
SIMV, SIMVASB (SIMV/PS).  
CPAP, CPAPASB (CPAP/PS).  
BIPAP1), (PCV+) (optional), BIPAP1)ASB (PCV+/PS) (optional).  
Enhancements - NIV – Non Invasive Ventilation with optimized alarm system and automatic leakage compensation (optional).  
Auto Flow®- Automatic adaptation of the inspiratory flow in volume orientated ventilation modes (optional).  
LPO - Low Pressure Oxygen. Independent oxygen supply, e.g. with an O2 concentrator (optional).  
Graphic screen - Advanced ventilation monitoring (optional).

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003