

Product Name :
Drager Savina

Product Code :
NLABMEDICAL-M3390010

Description :

Drager Savina

Technical Specification :

Ventilation modes - IPPV (CMV), IPPV Assist (CMV Assist).

SIMV, SIMVASB (SIMV/PS).

CPAP, CPAPASB (CPAP/PS).

BIPAP1), (PCV+) (optional), BIPAP1)ASB (PCV+/PS) (optional).

Enhancements - NIV – Non Invasive Ventilation with optimized alarm system and automatic leakage compensation (optional).

Auto Flow®- Automatic adaptation of the inspiratory flow in volume orientated ventilation modes (optional).

LPO - Low Pressure Oxygen. Independent oxygen supply, e.g. with an O2 concentrator (optional).

Graphic screen - Advanced ventilation monitoring (optional).

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003