

Product Name :
Draining Rack with Drip Pan

Product Code :
NLAB-EDUSCHOOLEQ207004



Description :

Draining Rack with Drip Pan

Technical Specification :

Depth: 200 mm

Width: 560 mm

Height: 680

White plastic-coated wire rack

37 rods with plastic caps

11runners

Plastic drip pan

Weight: 1,9 kg

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003