

Product Name :
Dressing Trolley Deluxe

Product Code :
NLABMEDICAL-M2810005



Description :

Dressing Trolley Deluxe

Technical Specification :

Frame Work of Dressing Trolley is made of a heavyweight Round tube.

It has stainless steel shelves with railing.

Dressing Trolley Deluxe is mounted on 4 castors of dia 10 cms.

It comes with a stainless steel bowl & bucket.

Overall Product Size : 90 L x 45 W x 80 H cms.

Knock Down Construction to save freight.

We are Manufacturers, Exporters, and Suppliers of Dressing Trolley Deluxe at tender and export competitive prices from India.

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments