

Product Name :
Dressing Trolley Standard

Product Code :
NLABMEDICAL-M2520001



Description :

Dressing Trolley Standard

Technical Specification :

Frame Work of Dressing Trolley is made of a Round tube.

It has stainless steel shelves with railing.

Dressing Trolley Standard is mounted on 4 castors of dia 5 cms.

It comes with a stainless steel bowl & bucket.

Overall Product Size : 90 L x 45 W x 80 H cms.

Knock Down Construction to save freight.

We are Manufacturers, Exporters, and Suppliers of Dressing Trolley Standard at tender and export competitive prices from India.

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments