

Product Name :
Dressing Trolley with Drawer

Product Code :
NLABMEDICAL-M2520003



Description :

Dressing Trolley with Drawer

Technical Specification :

Specification :-

Frame Work of Dressing trolley with drawer is made of heavy duty angles.

It has stainless steel shelves with up stands on four sides.

The Trolley is mounted on 10 cms castors.

It comes with stainless steel bowl size 32cm & Drawer.

The Dressing Trolley is fitted with one Drawer.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of Dressing Trolley with Drawer at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003