

**Product Name :**  
Drop Infusion Pump

**Product Code :**  
NLABMEDICAL-M3320001



**Description :**

Drop Infusion Pump India, Drop Infusion Pump manufacturer India, suppliers, Medical Drop Infusion Pump India, Hospital Drop Infusion Pump India

**Technical Specification :**

Offers a broad therapy spectrum with 0.5 to 5 liters O<sub>2</sub>/min.

Mobile dose monitor with its 20 meters long hose permits greater freedom of movement.

Very quite -only 40 dB.

Oxygen output more than 90%, outlet pressure 0.6 bar.

Alarms for mains failure/device failure and own O<sub>2</sub> concentration.

Small, light weight (only 20 kgs.), modular construction, snap closures, simplest filter change & low maintenance.

High performance and output even at 2000 meters above.

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## Naugralabequipments