

Product Name :
DROP-POINT APPARATUS

Product Code :
NLAB-ENGINEERINGLB16005



Description :

DROP-POINT APPARATUS

Technical Specification :

DESCRIPTION:-

This is used for determining drop point of greases as per IP 31,132 & ASTM-D-566. It is an indication of the critical temperature of the solid or semi solid to the liquid state. The apparatus consists of brass sleeve and case with metal cup and a glass boiling tube with cork fitted to a bath (Beaker) is provided a stirrer is also provided in the bath.

SERVICE REQUIRED:-

- Grease Sample.
- Suitable Table space to mount apparatus.
- 15 Amps. Single Phase Electrical Connection

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003