

**Product Name :**  
Dropping Bottles

**Product Code :**  
NLAB-EDUSCHOOLEQ131001



**Description :**

Dropping Bottles

**Technical Specification :**

Blown in Low Density Polyethylene, Dropping Bottles dispense small uniform drops of liquids when pressed lightly.

This Draining Rack provides an excellent option for post wash draining of tubes & bottles in the laboratory.

Dropping Bottles Capacity 60ml.

Dropping Bottles Capacity 125ml

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003