

Product Name :
Dropping Bottles

Product Code :
NLAB-EDUSCHOOLEQ131001



Description :

Dropping Bottles

Technical Specification :

Blown in Low Density Polyethylene, Dropping Bottles dispense small uniform drops of liquids when pressed lightly.

This Draining Rack provides an excellent option for post wash draining of tubes & bottles in the laboratory.

Dropping Bottles Capacity 60ml.

Dropping Bottles Capacity 125ml

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003