

Product Name :
Dry Cleaning Apparatus Testing Lab

Product Code :
NLAB-TESTINGLABM88001



Description :

Dry Cleaning Apparatus Testing Lab

Technical Specification :

Dry Cleaning Apparatus - Dry cleaning Apparatus is used for preparing fabric samples for flammability testing. This apparatus can also be used for testing color fastness and dimensional changes during dry cleaning and during laundering. & The Dry Cleaning Apparatus consists of a motor driven corrosion resistant stainless steel tumbler, which rotates around an axis 50°; inclined to tumbler axis at approximately 45 rpm. The tumbler is also provided with a stainless steel cover with a rubber seal and a quick acting positive clamp. The tumbler is rotated with the help of an electric motor through pulley-belt, reduction gearbox, and chain and sprocket drive

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003