

Product Name :
Drying Oven

Product Code :
NLAB-VOCATIONALTVET17014



Description :

Drying Oven

Technical Specification :

Capacity: 225 liters

Temp. range: +50-300°C

Time Locating Range: 1-9999 min

Temp. Fluctuation: ± 0.75 0C

Temperature Resolution: ± 0.10 C

Rated power heater elements: 2.5 KW

Shelves: 3

Shelf Positions: 4

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003