

<b>Product Name :</b> Dual Power Binocular Stereo Microscope	<b>Product Code :</b> NLAB-EDUSCHOOLEQ173098
<b>Description :</b> Dual Power Binocular Stereo Microscope	
<b>Technical Specification :</b> Durable and precise stereo microscope.  The fully coated optical system provides sharp and clear imaging with good flatness and contrast.  System Working Distance: 5.75-11.5mm  Standard Objective: 2/4X Objective  System Optical Magnification: 20X/40X  Optical System: Finite  This stereo microscope offers ergonomic features such as synchronized bi-lateral inter-pupillary distance adjustments and built in 3-way illumination.  Expandable System Working Distance: 120mm  Standard Eyepiece: 10X Eyepiece  Total Magnification: 20X/40X	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003