

Product Name :
Ductilomètre

Product Code :
NLAB-EDUCATIONTVET11006



Description :

Ductilomètre

Technical Specification :

Determination of the ductility of bitumens by measuring the elongation at the point of rupture of a bitumen test piece

Temperature: 25 ° C +/- 0.5 ° C

Maximum stroke: 1500 mm

Automatic feed: 50 mm / min

Equipped with a refrigeration unit for testing with water between 5 and 25 ° C

Can receive 3 test pieces simultaneously

Delivered without mold or base plate

Briquette mold EN 13398

Briquette mold EN 13589

Base plate

Ball and ring device

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003