

Product Name :
Dunagan Test Set

Product Code :
NLAB-VOCATIONALTVET17008

Description :

Dunagan Test Set

Technical Specification :

The balance of 5 kg capacity.

Containers: Immersion and soaking buckets with handles; water overflow container.

Density Basket: with dimensions 203 mm diameter x 203 mm deep, constructed of reinforced no. 8 wire mesh.

Weights: 2,000 g brass weight set and 1 kg slotted types.

Scoops: Set of three for weighing.

Sieves: 203 mm (8") diam.; one no.4 mesh x 4" (102 mm) h.; one no.100 mesh, standard height with no.100 mesh.

Case: Galvanized steel construction with cover and carrying handles.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003