

<b>Product Name :</b> Dynamo Bicycle	<b>Product Code :</b> NLAB-EDUSCHOOLEQ1730199
<b>Description :</b>	
Dynamo Bicycle	
<b>Technical Specification :</b>	
Two 4 mm socket terminals are provided for the dynamo output.	
Connected in parallel with an lamp holder and 2.5 V bulb.	
A cycle dynamo mounted on a wooden base 200 x 115 x 12 mm and provided with two sets of gears driven by a hand crank for fast or slow rotation.	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003