

Product Name : Dynamo Bicycle	Product Code : NLAB-EDUSCHOOLEQ1730199
---	--

--

Description :

Dynamo Bicycle

Technical Specification :

Two 4 mm socket terminals are provided for the dynamo output.

Connected in parallel with an lamp holder and 2.5 V bulb.

A cycle dynamo mounted on a wooden base 200 x 115 x 12 mm and provided with two sets of gears driven by a hand crank for fast or slow rotation.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003