

Product Name : Eccentric Presses	Product Code : NLAB-MACHINELABWK26004
--	---



Description :
Eccentric Presses
Technical Specification :
<p>Eccentric Presses .</p> <p>Indigenous make small eccentric presses specially designed for small pressing jobs such as stamping, marking, numbering, etc., ideally suitable for all small scale industries available in 3 . sizes Nos. 1, 2 and 3 represent small and big hand driven presses with vee drive arrangements :</p> <p>Approx. tonnage pressure Stroke of ram Stroke of back Front to back Base to guide Approx. net weight kgs.</p>

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003