

Product Name :
Edge Folding Machine

Product Code :
NLAB-MACHINELABWK24005



Description :

Edge Folding Machine

Technical Specification :

Edge Folding Machine .

oll type edge folding machines for edging the tin plate body blanks so form hooks for side seams prior to body forming or simple flanging with moderately sharp bends, complete with an adjustable gauge.

Length of fold

Capacity max. swg.

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com