

Product Name :
Educational Lab Bell Jar

Product Code :
NLAB-EDUSCHOOLEQ34016



Description :

Educational Lab Bell Jar


Technical Specification :

A bell jar is a piece of laboratory equipment used to contain vacuum.

It is similar in shape to a bell, and can be manufactured from a variety of materials ranging from glass to different types of metals.

The base of the jar is equally heavy and flattened.

Glass, with ground glass stopper at top, ground rim.



226

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003