

Product Name :
Educational Lab Bicycle Dynamo

Product Code :
NLAB-EDUSCHOOLEQ34018



Description :

Educational Lab Bicycle Dynamo

Technical Specification :

A bicycle dynamo is a model that converts mechanical motion into electric motion through a rotating wheel and magnet.

The dynamo has a perpetual magnet, twisted with coils of cloistered wire.

When the magnet rotates, it creates a magnetic field, which in turn creates electricity in the wire.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003