

Product Name :
Educational Lab Cryo Vials

Product Code :
NLAB-EDUSCHOOLEQ34058



Description :

Educational Lab Cryo Vials

Technical Specification :

Cryo vials are normally used for storage of Biological material, human & animal cells.

These sterilized tubes are provided with white, prominent graduations and a long marking area for sample identification.

The cap of these vials is provided with a seal ring inside to make it completely leakproof.

The longitudinal grooves in the base give a non-twisting feature to these vials when placed in a cryo rack.

Material : Polypropylene

Specification:

Cryo Vail 1.0 ml

Cryo Vail 1.8 ml

Cryo Vail 4.5 ml

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003