

Product Name :
Educational Lab Desoldering Pump

Product Code :
NLAB-EDUSCHOOLEQ34063



Description :

Educational Lab Desoldering Pump

Technical Specification :

As the name suggests, desoldering pump is a suction unit that sucks the molten solder when the button is triggered and can be removed later from the device.

It is made from premium materials to endow a better performance and functioning to the pump.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003