

**Product Name :**  
Educational Lab Digital Altimeter

**Product Code :**  
NLAB-EDUSCHOOLEQ34066



**Description :**

Educational Lab Digital Altimeter

**Technical Specification :**

Compact.

Accurate results.

Wide measurement range.

Low priced model for personal use.

Ideal for hiking, camping, mountaineering etc.

6-in-1 multi-purpose model with height, pressure, direction, temperature, weather forecast, time & date.

Specifications:

Altimeter range: -2296 ft to 29500 ft (-700m to 9000m)

Resolution: 1m/3ft

Measurement: sea level altitude/climbing altitude

Unit: M and Ft switchable

Barometer range: 300hpa to 1100hpa (8.8 to 32.6inHg)

Resolution: 0.1hpa

Unit: Millibar and INhg switchable

Thermometer range: -20 to 60C (-4 to 140F)

Resolution: 1C/1F

Unit: Celsius and Fahrenheit switchable

Compass range: 360 degrees with 8 directions

Resolution: 1 degree

Weather Forecast: 12-24 hours; sunny-slightly sunny-cloudy-rain

Time & Date: 24H or 12 H

Back-light: 5 seconds super bright LED

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003