

**Product Name :**  
Educational Lab Digital Conductivity Meter

**Product Code :**  
NLAB-EDUSCHOOLEQ34067



**Description :**

Educational Lab Digital Conductivity Meter

**Technical Specification :**

Display 3½ LED.

Range 0 to 1000 Mhos/cm 5 Range.

Accuracy  $\pm 1\%$  FS  $\pm 1$  Digit.

Power Requirement 230V AC.

Spare Conductivity Cell 1./0.5/0.1.

Complete with Operation Manual & Dust Cover.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

