

Product Name :
Educational Lab Digital Penetrometer

Product Code :
NLAB-EDUSCHOOLEQ34070



Description :

Educational Lab Digital Penetrometer

Technical Specification :

Digital Penetrometer is used for penetration test on different materials.

It contains a head mounted on a perpendicular rod which can be attuned for height.

Penetrative Strength of the material is determined with help of this apparatus.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003