

Product Name :
Educational Lab Draining Rack

Product Code :
NLAB-EDUSCHOOLEQ34076



Description :

Educational Lab Draining Rack

Technical Specification :

Draining rack provides an excellent option for post wash draining of tubes & bottles in the laboratory.

The back plate has 20 holes into which pegs are tightened with the help of screw nuts.

Bottom of the rack is provided with a tray which collects all the waste water drained down which can further be collected in a container via drain tube.

This draining rack does not rot, stain or rust.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003