

Product Name :
Educational Lab Elastic Material Kit

Product Code :
NLAB-EDUSCHOOLEQ34082



Description :

Educational Lab Elastic Material Kit

Technical Specification :

In this kit, various materials such as springs, rubber blocks, wire, rubber cord, tube etc., are provided for class experiments to explore 'elasticity' by performing experiments on compression, torsion, extension and elastic limits.

The standard pack is the one-quarter class kit, i.e. sufficient for eight students and contains:

25 expendable steel springs.

2 wide steel springs 114 x 75 mm diameter.

2 soft latex foam blocks.

1 50g reels of 32 s.w.g. Bare copper wire.

2 50g reels of 26 s.w.g. Bare copper wire.

1 length valve rubber tubing.

4 elastic cords with eyelets.

4 soft erasers for twisting.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003