

Product Name :
Educational Lab Galvanometer, DC

Product Code :
NLAB-EDUSCHOOLEQ34097



Description :

Educational Lab Galvanometer, DC

Technical Specification :

Galvanometer, DC is a device used for measuring small electric currents.

When a current passes through the device, the coil feels a torque proportional to the charge.

The device is fitted with a meter in which the current is shown through a moving needle.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003