

**Product Name :**  
Educational Lab Galvanometer, DC

**Product Code :**  
NLAB-EDUSCHOOLEQ34097



**Description :**

Educational Lab Galvanometer, DC

**Technical Specification :**

Galvanometer, DC is a device used for measuring small electric currents.

When a current passes through the device, the coil feels a torque proportional to the charge.

The device is fitted with a meter in which the current is shown through a moving needle.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003