

Product Name : Educational Lab Ganong Potometer	Product Code : NLAB-EDUSCHOOLEQ34099
Description : Educational Lab Ganong Potometer	
Technical Specification : Ganong Potometer holds a graduated narrow bore glass tube attached with a vessel to carry a plant shoot at one end and funnel with glass stopcock on other end. It is used for measuring the rate of water uptake of a leafy plant shoot. It comes complete with wooden stand.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003