

**Product Name :**  
Educational Lab Ganong Potometer

**Product Code :**  
NLAB-EDUSCHOOLEQ34099



**Description :**

Educational Lab Ganong Potometer

**Technical Specification :**

Ganong Potometer holds a graduated narrow bore glass tube attached with a vessel to carry a plant shoot at one end and funnel with glass stopcock on other end.

It is used for measuring the rate of water uptake of a leafy plant shoot.

It comes complete with wooden stand.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003